

Mercury® Pylon

SHOCK ABSORPTION and AXIAL ROTATION

The new Mercury® Pylon from Endolite reduces vertical shock and permits controlled axial rotation between the socket and the foot. A larger bearing surface gives a smoother, more consistent and durable operation.

- A 35 mm pylon with adapters available for use with a variety of 35 and 30 mm systems.
- 30 degrees of rotation in either direction from the neutral position.
- 1/5 — 5/16 inch vertical travel during normal walking; 1/2 inch maximum.
- The latest torsion rod does not require adjustment.
- Can be used with a variety of feet.
- Reduced weight, longer bearing surface, silent operation and smoother function.

Innovation in stride

smooth

Circle # 58 on Readers' Service Card

ENDOLITE
Innovation in stride.™

mercury®