

# Stance - Control Orthoses™

Designed for people who have knee instability due to quadriceps weakness or absence. Results in a more normal and energy-efficient gait than traditional locked knee joints.



Polio Patient wearing Bi-lateral SCOKJ's®

**Horton Technology Inc.**

*Innovations In Orthotics and Prosthetics*

Prevents knee collapse regardless of flexion

Stance control can be triggered by weightbearing or joint motion, according to patient needs

Provides three operating modes: Automatic Stance Control, Free Motion, and Manual Lock

Can be used with bilateral knee paralysis or paresis

[www.stancecontrol.com](http://www.stancecontrol.com) 1-866-663-3970