



M-**BRACE**™

soft wear for the body



BRACE YOURSELF

BACK BRACES #582 & #584

This LumbLock Lumbar-Sacral Brace acts as an immobilizing device for the lumbar-sacral region and its' primary function is to reduce pain caused by lumbago, sciatica and ischialgia.

Made with semi-rigid cotton fabric, it enjoys a structure made with eight selfmodeling and flexible steel stays on which a number of rigid straps slide. When the locking system is pulled and secured at the front, the inside straps stiffen the entire structure, thus effectively immobilizing the lumbar-sacral region.

Practical and easy to wear, this LumbLock brace perfectly adapts to whatever size and body type, thus ensuring maximum support and comfort at the desired tension.

Approved for L-0637, L-0631, L-0627

Use these back braces for the following indications:

- Postural pathologies: Lumbago, sciatica, discopathy, hypo-lordosis. (582)
- Supports from sacrum to T-9 level. (584)
- Post-op: Following stabilizing surgery, following lumbar disk surgery. (584)
- Degenerative pathologies: Osteoporosis, lumbar- sacral vertebrae complaint, protrusions, arthritis. (584)
- Degenerative pathologies: Osteoporosis, vertebrae complaints, protrusions, arthrosis. (582)
- Muscular pathologies: Muscular contractures.
- Practical and easy to fit, it adapts perfectly to your body.
- Breathable and non allergenic



Use this knee brace for the following indications:


- Mild to moderate medial compartment OA without instability.
- Meniscus procedures.
- Chondral defect procedures.
- Moderate lateral compartment OA without instability; it is possible to reduce the load on the lateral compartment by using a "right" brace on a left leg or a "left" brace on a right leg.
- Available in left or right.

KNEED SUPPORT?

OA KNEE BRACE WITH RANGE OF MOTION

In your quest to live an active lifestyle, the #46 knee brace brings your comfort to a new level with its elegant and light design. Available for the right or left knee, this brace alleviates knee pain by reducing the weight on your medial or lateral compartment by using the three point pressure system. It is composed of two aluminum reinforced semi-rigid straps that surround the distal region (above the knee) and the proximal region (below the knee). These straps are to be custom moulded to allow for an optimal torsion support. Additionally, these straps are bound together by a patented articulated aluminum support that is easily adjustable for flexion and extension and accessible within a zipped pocket.

KNEE #46



Approved for L-1843

Contact us today for more information.



www.m-brace.com

1 866-374-4050