

# Imagine the possibilities!

The Motionfoot® truly improves the comfort and naturalness of walking and using an artificial foot. These features include:

- A hydraulic ankle with up to 50 degrees of motion, allowing natural-like plantar flexion (pointing the toe downward) and dorsi flexion (pointing the toe upward)
- The wearer can adjust the resistance of the ankle. Two independent hydraulic settings allow you to match the ankle to your activities.
- The wearer can also manually lock the ankle at 90 degrees (neutral standing position). This allows you to prevent unwanted plantar flexion for specific activities.
- The carbon fiber Sole Plate can also be adjusted (by your prosthetist) to match your weight and activity level.

**Motionfoot**®

*MotionControl*  
KEEPING LIFE IN MOTION

115 North Wright Brothers Drive • Salt Lake City, Utah 84116  
801-326-3434 • 888-696-2767 (toll free) • 801-978-0848 (fax)  
[www.UtahArm.com](http://www.UtahArm.com)