

# As essential as morning coffee



## The Otto Bock Harmony® Volume Management System

Starting your patient's day with the uniquely effective Harmony® Volume Management System is as comforting as a morning cup of coffee.

By adding the Harmony System to your patient's prosthesis, you apply vacuum to maintain the ideal socket environment. The results are much more secure suspension and all the benefits of improved connectivity. The Harmony System also delivers volume management – a great decrease in the variations driven by fluid changes in the residual limb – meaning that your patient gets improved fit without the hassle of adding and removing socks.

To add Harmony to your patient's mornings, simply call us at 800.328.4058 (US) or 800.665.3327 (Canada) or visit [www.ottobockus.com](http://www.ottobockus.com) for more information and course listings.