

## Different fits for different folks.



The Harmony® Volume Management System now comes in three styles: The new Harmony DP is designed exclusively for use with the Advantage DP pylon foot, the Harmony P2 for patients weighing between 110–220 lbs and the Harmony HD for those weighing between 220 lbs and 330 lbs.

The Harmony system uses a pump, liner and sleeve to achieve elevated vacuum in an airtight environment. Add a professionally designed socket, and your patients benefit from incredible linkage, improved proprioception and better tissue health.

To learn more about delivering the best fit to all your patients, contact Otto Bock at **800.328.4058** (US) or **800.665.3327** (Canada).