



HOPe1 Hip Orthosis Pediatric 1

Indications:

Night bracing post multi-level Botox in spastic cerebral palsy (quadriplegia, diplegia), or where abduction of hip is needed.

Features and functions (benefits):

- 1) Standard Ultraflex platform (joint) and power unit mounts to medial mid calf KO or KAFO sections with swivel brackets:
 - Provides 7.5 degree increment don/doff locks
 - 7.5 degree increment adduction or abduction stop (please specify)
 - Adjustable dynamic tension for abduction assist
 - Adjustable internal/external rotation position
 - Components gently abduct to stable full maximum abduction
 - Lifetime limited component warranty on original patient user
- 2) Full "gutter" solid KO or KAFO sections:
 - Minimize point pressure
 - Simple two strap patellar counter force strapping keeps knee in full extension (as cast)
 - Can be used with abduction bars/component removed as individual KO's or KAFO's
 - Can be used with articulated designs (knee/ankle)



Easier to don providing better compliance than traditional "A" Frames!



Quick release components allow more flexibility of use

UltraflexSystemsInc.

237 South Street • Suite 200
Pottstown, Pennsylvania 19464

800-220-6670

www.ultraflexsystems.com

Copyright 2005 Ultraflex Systems, Inc.

UltraflexEurope

distributed by
DIRAME ORTHO nv/sa
Assesteenweg 27-29
B-1740 Ternat (Brussels) - Belgium
+32/(0)2/582.82.50
www.dirame.be