

genu range of motion

② $2 \frac{2}{3} \div \frac{2}{3}$

CONFIDENCE

to stand at every angle.

Increase →



Circle # 165 on Readers' Service Card



STEP LOCK

Orthotic Knee Joints

- > Ten locking positions throughout knee flexion
- > Helps to increase the knee range of motion
- > Provides stability for standing before reaching full extension
- > Designed for knee flexion contracture management
- > Unlocks to enable sitting and reclining
- > Effectively blocks flexion while allowing free extension
- > Unparalleled strength plus a nickel plated finish for extra durability

OTS Corp. • 1.800.221.4769

Fillauer LLC • 1.800.251.6398

Centri® • +46 8 505 332 00

Hosmer • 1.800.827.0070